Thermo Scientific Nunc FluoroNunc Plates and Modules: A Solid Phase for Fluorescent Immuno Assays

Key Words

Thermo Scientific[™] Nunc[™] FluoroNunc[™], Thermo Scientific[™] Nunc[™] PolySorp[™], Thermo Scientific[™] Nunc[™] MaxiSorp[™], solid phase,

Goal

The goal of this Tech Note is to describe the use of white Nunc FluoroNunc C96 plates and white C8 modules in a fluorescence immuno assay.

This Tech Note describes the use of white Nunc FluoroNunc C96 plates and white C8 modules in a fluorescence immuno assay.

The use of fluorescence techniques in solid phase diagnostic assays is attracting increasing attention due to the higher sensitivity which can potentially be obtained. Sensitivity, however, is not only a question of the right detection system, but is also dependent on the performance of the solid phase used.

The Nunc FluoroNunc plates and modules have been developed in a transparent version for use in time resolved fluorescence assays where the plate background is minimized. The black version minimizes light scattering, while the white version is used for high sensitivity in traditional fluorescence assays (Tech Note No. 6).

All the Nunc FluoroNunc products can be obtained with Nunc PolySorp or Nunc MaxiSorp surface which provides the binding properties needed to make the most sensitive and reliable assay.

To demonstrate the use of white Nunc FluoroNunc plates and modules we have designed a sandwich assay for detection of IgG.

Materials and Method

Solid phase: Nunc FluoroNunc Plate, White C96 Nunc MaxiSorp.

Antibodies: Pig anti-rabbit, Rabbit anti-sheep, Pig anti-rabbit alkaline phosphatase conjugate.

Buffers: Carbonate buffer 0.05M pH 9.6. Phosphate buffered saline 0.15M pH 7.2 (PBS). Assay buffer: PBS with 0.05% Tween 20. Wash buffer: PBS with 0.2 M NaCl and 0.05% Triton X100.



Substrate: 0.1 mm 4-methyl-umbelliferylphosphate in diethanolamine buffer pH 9.8.

Miscellaneous: Bovine serum albumin.

Stop solution: 3 M K₂HPO₄.

150 µL carbonate buffer containing 5 µg/mL PaR antibody was added to each well. The wells were sealed with adhesive foil. The plates were coated overnight at room temperature. Then the plates were washed once with demineralized water. Blocking solution (400 µL PBS with 0.5% BSA) was added to all the wells and the plates incubated for 10 minutes at room temperature. The wells were completely aspirated, and a 1:2 dilution series of RaS, 150 µL/well, was performed in the plates (starting in column 2 at a concentration of 420 ng/mL). Assay buffer was added to column 1. The plates were incubated for 2 hours at room temperature, then washed 5 times with wash buffer using Thermo Scientific[™] Nunc[™] Immuno Washer.



Conjugate in assay buffer, 150 μ L/well, was added to all wells in dilution 1:500 or 1:5000, and the plates incubated for 2 hours at room temperature. After the plates had been washed 5 times with wash buffer the substrate reaction was initiated by adding 150 μ L/well of substrate. The reaction was stopped after 15, 30 or 45 minutes, by adding stop solution.

The fluorescent results were read on a fluorescence reader and the data collected on computer using in house software.

Results

The results of this test are shown in Fig. 1 and graphically in Fig. 2.

As can be seen the best sensitivity was obtained using conjugate in a dilution of 1:500 and a substrate reaction time of 30 minutes. No improvement in sensitivity was observed when the substrate reaction time was increased to 45 minutes. Data analysis, using background (Bg) signal plus 2 SD as cut off value showed that a detection of 20-40 pg/mL RaS antibody could be clearly obtained in this test.

Fig. 1

Results of fluorescence assay

15 min; 1:500

ng/mL lgG	0.01	0.02	0.04	0.08	0.16	0.31	0.625	1.25	Bg
Mean	26.77	25.41	38.98	54.39	82.12	148.2	274.3	551.3	8.91
SD	11.11	5.253	13.45	14.24	7.941	6.554	16.32	28.32	4.77
2 SD	22.22	10.506	26.9	28.48	15.882	13.108	32.64	56.64	9.548
Mean + 2 SD	48.99	35.916	65.88	82.87	98.002	161.308	306.94	607.94	18.459
Mean - 2 SD	4.55	14.904	12.08	25.91	66.238	135.092	241.66	494.66	

30 min; 1:500


ng/mL lgG	0.01	0.02	0.04	0.08	0.16	0.31	0.625	1.25	Bg
Mean	41.03	53.07	76.4	116.7	179.2	315.6	614.2	1145	21.21
SD	10.44	9.865	12.42	29.56	9.282	22.68	20.56	59.57	5.594
2 SD	20.88	19.73	24.84	59.12	18.564	45.36	41.12	119.14	11.188
Mean + 2 SD	61.91	72.8	101.24	175.82	197.764	360.96	655.32	1264.14	32.398
Mean - 2 SD	20.15	33.34	51.56	57.58	160.636	270.24	573.08	1025.86	

45 min; 1:500

ng/mL lgG	0.01	0.02	0.04	0.08	0.16	0.31	0.625	1.25	Bg
Mean	66.38	84.29	131	194.5	298.6	542	1018	1949	56.28
SD	18.35	21.22	38.62	50.23	17.88	40.8	39.1	75.32	16.54
2 SD	36.7	42.44	77.24	100.46	35.76	81.6	78.2	150.64	33.08
Mean + 2 SD	103.08	126.73	208.24	294.96	334.36	623.6	1096.2	2099.64	89.36
Mean - 2 SD	29.68	41.85	53.76	94.04	262.84	460.4	939.8	1798.36	

Fig. 2

Graphical representation



Discussion

The successful performance of solid phase fluorescent technique is dependent on the solid phase used just as much as on a reliable detection system. The result of this test has shown that a sensitive fluorescent assay can be obtained on Nunc FluoroNunc white plates and modules. Furthermore, these products have the same quality of surface performance as all Nunc Immuno products. Existing chromogenic tests designed on Nunc PolySorp or Nunc MaxiSorp surfaces can therefore easily be transferred to fluorescence tests by using Nunc FluoroNunc plates or modules and simply changing the substrate to a fluorescent substrate offering an increased assay sensitivity.

thermoscientific.com/oemdiagnostics

© 2014 Thermo Fisher Scientific Inc. All rights reserved. "Tween" is a registered trademark of Uniqema Americas; and "Triton" is a registered trademark of Dow Chemical Company. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

ANZ: Australia: 1300 735 292, New Zealand: 0800 933 966; Asia: China Toll-free: 800-810-5118 or 400-650-5118; India: +91 22 6716 2200, India Toll-free: 1 800 22 8374; Japan: +81-3-5826-1616; Other Asian countries: 65 68729717 Europe: Austria: +43 1 801 40 0; Belgium: +32 2 482 30 30; Denmark: +45 4631 2000; France: +33 2 2803 2180; Germany: +49 6184 90 6000, Germany Toll-free: 0800 1-536 376; Italy: +39 02 95059 554; Netherlands: +31 76 571 4440; Nordic/Baltic countries: +358 9 329 10200; Russia/CIS: +7 (812) 703 42 15; Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 22; UK/Ireland: +44 870 609 9203 North America: USA/Canada +1 585 586 8800; USA Toll-free: 800 625 4327 South America: USA/Sales support: +1 585 899 7198 Countries not listed: +49 6184 90 6000 or +33 2 2803 2000