

ROTO SILENTA 630 RS Floorstanding Centrifuge cooled



ROTO SILENTA 630 RS

The ROTO SILENTA 630 RS offers outstanding performance in fast, efficient centrifuging of blood bags and large sample volumes. It accommodates a maximum capacity of 6 x 2,000 ml or 12 x 750 ml blood bag systems.

Excellent stability even at high speeds makes the ROTO SILENTA 630 RS a rugged and reliable centrifuging partner. When in use, the cooled large-volume centrifuge rests on high-quality vibration dampers.

To ensure mobility, it is supplied with castors onto which the ROTO SILENTA 630 RS can be lowered, as required.

Optimal running performance guarantees best separation results. The low noise level makes agreeable working conditions. The extensive range of carriers and adapters for blood collection tubes, standard tubes, bottles and various blood bag systems can be centrifuged in one rotor – changing rotors is not necessary.

The Hettich data report system HettInfo as well as customised solutions enable GMP performance and meet the high standards of user safety and product quality in blood banks and transfusion centres.

PERFORMANCE

- High RCF without windshield
 - · up to 6,520 x g at 4,500 min⁻¹.
- High capacity
 - · max. 6 x 2,000 ml or 12 x 750 ml blood bag systems
- Efficient cooling
 - \cdot controllable from 20 °C to + 40 °C (with model 5005-50 from 20 °C to + 90 °C).
- Precooling is effected by programming one of the 89 programmable memories.

TECHNOLOGY

- Frequency-controlled drive, virtually maintenance-free
- Powered lid locking
- Ergonomic rail-grip for easy lid closing
- Large viewing port in the lid for optical speed control with a stroboscope

SAFETY

- Imbalance switch-off
- Rotor recognition
- Key-operated switch for locking and protecting programs
- Lid locking and holding
- Lid dropping protection
- Motor and chamber overheating protection
- Torsion-free steel chassis
- Smooth, easy to clean surfaces



ROTO SILENTA 630 RS with bar code unit (option), dimensions (H x W x D): 973 x 813 x 1,015 mm (without bar code unit), 973 x 888 x 1,015 mm (with bar code unit). On request, the scanner can also be mounted on the left panel.

Cat.-No. 5005

CONTROL PANEL

ROTO SILENTA 630 RS with S control panel to satisfy the highest standards

Numerous functions with just a few user operations

Whether entering single parameters or chaining complete running programs, the clearly arranged display with its selector keys and adjusting knob provides for extreme user-friendliness. The controls and displays of the ROTO SILENTA 630 RS are arranged according to ergonomic aspects. The large display is clearly arranged and well readable. During centrifugation the

For special routine or scientific applications the ROTO SILENTA 630 RS can be fitted with the following optional equipment:

actual values of all parameters are indicated.

Centrifuging data report system HettInfo details on next page

Visual indication when the rotor is stationary

for additional information on the completion of a centrifuging cycle.

Program interlocking

to combine several centrifugation runs.

Second independent tacho system

for calibrating the speed measuring instrument as required by quality management systems.



ROTO SILENTA 630 RS control panel: Entry of the parameters is exact, fast and easy with selector keys on the foil keypad and the adjusting knob.

KEYPAD

START Starts centrifugation run, short-time centrifugation.

Pre-cooling can be assigned to one program location.

STOP Stops centrifugation manually.

PROG Selects the program.

89 programmable memories are available.

STO Stores entries, changes and programmes.

RCL Calls up the selected program.

JRCF The integral over the RCF indicates the

accumulated RCF acting on the centrifuged

material up to that particular time.

PARAMETER INPUT

t

Entry of the temperature from –20 °C to +40 °C in increments of 1 °C.

Entry of the centrifuging radius in mm.

Entry of the run-up time in ramps 1 – 9 or in min: sec.

Entry of the run-down time in ramps R 1-9

and B 1-9 or in min: sec.

Also unbraked run-down (setting 0),

or a brake force cut-off speed can be selected.

RCF Entry of the RCF in increments of 1.

n Entry of the speed in increments of 10.

Entry of the centrifugation time (max. 999 min: 59 sec)

or continuous operation.

Cool solutions under any conditions

The ROTO SILENTA 630 RS adapts to the most diverse local conditions and requirements. Depending on fittings, it can be connected to

an external cooling an external cooling unit medium circuit cooling unit integrated

green = cooling medium (R 404A), blue/red = coolant

Of course, a stand-alone model is also available: an air-cooled standard version supplied with an integrated cooling unit.


The housing of all versions can be insulated for use under clean-room conditions.

Additionally, the pedestal can be covered and sealed to the floor. External connections can also be mounted through the bottom of the centrifuge.

In addition to the above mentioned centrifuge models, we also manufacture special versions and accessories tailored to customer requirements.

Transparent and reliable processes thanks to HettInfo

With the data documentation system HettInfo, up to 29 centrifuges can be controlled by a single PC. A light conducting cable connects the centrifuge to a BUS system. A PC with a RS 232 interface records the selected data of the centrifuge.



Hettlnfo registers all centrifugation data. In addition to the date, the end of the centrifugation run, the program number and the temperature, further parameters can be selected or deselected by the user. Detailed information is available in our Hettlnfo brochure. Please ask for it.



Packing room with external bar code unit



Barcode-Kit ROTO SILENTA 630 RS



Centrifugation programs and the operator code can also be scanned-in by barcode.

Data log

MODEL VARIANTS AND ROTORS

Model variants ROTO SILENTA 630 RS

Variant	Power supply*)	Frequency	Cat. No.
With connection for nitrogen (N ₂) flushing	400 V 3~+N	50-60 Hz	5005-20
Heating/cooling version up to + 90 °C	400 V 3~+N	50-60 Hz	5005-50
Brine cooling	400 V 3~+N	50-60 Hz	5005-60
GMP, internal cooling unit with water-cooled condenser	400 V 3~+N	50-60 Hz	5005-80
GMP, external cooling unit	400 V 3~+N	50-60 Hz	5005-90

Rotors

With the ROTO SILENTA 630 RS there is a choice between a 4-place and a 6-place swing-out rotor. Both rotors can hold blood bags, standard tubes, blood collection tubes and bottles with a capacity of up to 2,000 ml.

The ROTO SILENTA 630 RS generates a maximum RCF of 6,520 with the 6-place rotor 4176 and a maximum RCF of 5,705 with the 4-place rotor 4174 – in each case **without windshield**. This not only simplifies handling, it also enables optimal temperature distribution and precise temperature measurement.

Lower temperatures can be achieved with 4-place rotor 4174 than with rotor 4176 thanks to its lower air friction values.

The values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background in the following tables of accessories.

Swing-out rotor, 6-place, ≮90°



 $n = 4,500 \text{ min}^{-1},$ max. RCF 6.520



illustrated with carriers 4595-C and bottles 0550 for 2,000 ml

n = 3,500 min⁻¹, max. RCF 3,848



illustrated with carriers 4522-A and two each inserts 4245-A, 4232 and 4223

 $n = 4,500 \text{ min}^{-1},$ max. RCF 5,999

Cat. No. (without carriers) 4176

^{*)} Other voltages on request.

Swing-out rotor, 6-place

In the following tables, the values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background.



illustrated with carriers 4591-A and inserts 4592-B

Cat. No. (without carriers) 4176

capacity in ml	500	500	500	750		
holding system for blood bags	-	-	-	-		
blood bags	quadruplicate	quadruplicate	quadruplicate	single		
Cat. No.	•	-	-			
carrier Cat. No. 4524-A						
Cat. No.	4525-A	4529 ⁴⁾ -AO,-AM,-AU	459	4592-B		
boring Ø x L in mm		-				
blood bag systems per rotor		12 8				
max. RCF		6,498 5,683				
max. RCF ⁵⁾		1,000 1,000				
radius in mm		287 251				
run-up in sec		125 125				
run-down in sec, braked		197 197				
temperature in °C1)	_	+16 +10				

capacity in ml	450						
blood bags, quadru	plicate						
Cat. No.	-						
8							
carrier Cat. No. 4546-A							
Cat. No.	4559-A						
boring Ø x L in mm	-						
blood bag systems per rotor	12 8						
max. RCF	6,271 5,479						
radius in mm	277 242						
run-up in sec	125 I 125						
run-down in sec, braked	197 <mark>197</mark>						
temperature in °C1)	+13 0						



carrier 4524-A with insert 4525-A

insert 4529-A in 4529-AM version (hooks for suspending the blood bags are fixed in middle position)



carrier 4546-A with insert 4559-A

¹⁾ Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be reached by reducing the speed.

³⁾ As well as Macopharma and Kawasumi.

When using inserts 4529-AU, 4529-AU at lower speeds (e.g. for processing thrombocytes), the blood bags may be hung up to avoid accumulation of erythrocytes. Hooks at three different heights (4529-AO upper position, 4529-AM middle position, 4529-AU low position) allow to meet specific customer requirements. Suspended blood bags, however, may not be exposed to a RCF exceeding 1,000.

ROTORS AND ACCESSORIES

Swing-out rotor, 4-place



illustrated with carriers 4591-A and inserts 4592-B max. RCF 5,705

Cat. No. (without carriers) 4174

capacity in ml	450	500	750			
blood bags	triplicate	single				
Cat. No.	-		-			
Carrier Cat. No. 4591-A						
Cat. No.	4598	459	2-B			
boring Ø x L in mm		-				
blood bag systems per rotor		12 8				
max. RCF		6,498 <mark>5,705</mark>				
radius in mm		287 <mark>252</mark>				
run-up in sec		125 125				
run-down in sec, braked		197 <mark>197</mark>				
temperature in °C1)		+16 +10				

capacity in ml	max. 1.600 ⁵⁾	2.000
	-	180x150x100
Cat. No.	-	0550
Carrier Cat. No. 4595-C ⁶⁾		
Cat. No.	4596-A	-
(HxWxD) in mm	146 x 145 x 91	-
blood bag systems per rotor	6 l 4	6 I <mark>4</mark>
max. RCF	3,821 <mark>3,328</mark>	3,848 <mark>3,369</mark>
radius in mm	279 I <mark> 243</mark>	281 <mark>246</mark>
run-up in sec	125 I <mark>125</mark>	29 I <mark>95</mark>
run-down in sec, braked	197 <mark>197</mark>	41 131
temperature in °C1)	+16 I <mark>-8</mark>	+16 I <mark>-8</mark>



carrier 4591-A with insert 4592-B



carrier 4595-C with insert 4596-A

Dependent on the material to be centrifuged.
 Adapter for accommodating sample tubes and blood collection tubes in carrier 4595-C on request.

In the following tables, the values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background.

Swing-out rotors 4176 / 4174

capacity in ml	500	750	1.000		
blood bags	quadrupli- cate	single	single		
Cat. No.			-		
Carrier Cat No. 4523-A					
Cat. No.	451	6-A	-		
boring Ø x L in mm			-		
blood bag systems per rotor	6	4	6 4		
max. RCF	6,452	5,637	6,520 5,705		
radius in mm	285	249	288 252		
run-up in sec		125	125		
run-down in sec, braked		197	197		
temperature in °C1)		+9	+3		

capacity in ml	500	450				
blood bags	quadrupli- cate	triplicate				
Cat. No.		•				
lid Cat. No. 5621						
Carrier Cat No. 4547-B						
Cat. No.	4548					
boring Ø x L in mm						
blood bag systems per rotor	6	4				
max. RCF	6,316	5,501				
radius in mm	279	243				
run-up in sec	125	125				
run-down in sec, braked	197	197				
temperature in °C1)		-9				



carrier 4523-A with insert 4516-A



carrier 4547-B with insert 4548 and lid 5621

¹⁾ Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be reached by reducing the speed.

ACCESSORIES FOR CENTRIFUGING BLOOD BAG SYSTEMS



Quick and easy packing of blood bag systems with loading aid

Cat. No. 4509



HettLiner

Blood bags can be secured in adapters using HettLiners. After centrifuging, they can be removed by their loops without disturbing the sediment.

Cat. No. 4564



Balancing inserts

In case there are not enough blood bag systems to occupy every carrier of the rotor, empty carriers can be filled with balancing inserts. Taring weights supplied with the inserts may be used for fine balancing.

Cat. No.

4584-A (for insert 4559-A)

4587-A (for insert 4592-B)

4589-A (for insert 4516-A)



Set of compensatory weights

If blood bag carriers are not filled to the same weight, then compensatory weights can be used to make up the difference.

Cat. No. 4566

In the following tables, the values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background.



carrier 4579-A, illustrated with bucket 4255 with lid

Swing-out rotors 4176 / 4174

capacity in ml	5	7	9	15	25	50	100	2.6-3.4	2.7-3	
Ø x L in mm	12x75	12×100	14×100	17×100	24×100	34×100 44×100		13×65	11×66	
Cat. No.	0553 ²⁾	0578 ²⁾	0500 ²⁾	0518 ²⁾	0519 ²⁾	0521 ²⁾	0526 ²⁾	blood colle	d collection tubes	
bucket Cat. No. 4255	t. No. 4255									
Carrier Cat. No. 4579-A					+ 0726					
Cat. No.	44	33	44	34	4438	4439	4442	4435	4433	
boring Ø x L in mm	13:	x58	17.5	x53	26x72	36×79	45×78	13.5 x 58	13x58	
tubes per rotor	180	120	114	76	42 <mark>28</mark>	24 <mark>16</mark>	12 <mark> 8</mark>	126 <mark>84</mark>	180 120	
max. RCF ²⁾	5,750	4,935	5,886	5,094	5,615 <mark>4,845</mark>	5,705 4,890	5,683 4,867	5,750 4,935	5,750 4,935	
radius in mm	254	218	260	225	248 <mark>214</mark>	252 216	251 <mark>215</mark>	254 <mark>218</mark>	254 218	
run-up in sec			•		125	125				
run-down in sec, braked					197	197				
temperature in °C1)					+4	-11				

capacity in ml	4-5.5	4.5-5	4.9	9-10	10	1.6-5	4-7	4-7	8.5-10	15	50	
Ø x L in mm	15×75	11x92	13×90	16x92	15×102	13x75	13×100	16x75	16 x 100	17×120	29×115	
Cat. No.	blood collection / urine tubes										0513	
bucket Cat. No. 4255					7)	J						
Carrier Cat. No. 4579-A												
Cat. No.	4434	4433	4435	44	34	44	35	4434		4437	4441	
boring Ø x L in mm	17.5 x 53	13×58	13.5 x 58	17.5	x53	13.5	x58	17.5×53		17×88	30×87	
tubes per rotor	114 <mark>76</mark>	180 120	126 84	114	76	126	84	114	76	72 <mark>48</mark>	30 20	
max. RCF ²⁾	5,886 <mark>5,094</mark>	5,750	4,935	5,886	5,094	5,750	4,935	5,886	5,094	6,022 <mark>5,207</mark>	5,999 <mark>5,207</mark>	
radius in mm	260 <mark>225</mark>	254	218	260	225	254	218	260	225	266 <mark>230</mark>	265 230	
run-up in sec					125	1 125						
run-down in sec, braked					197	197						
temperature in °C1)					+4	ı -11						

10

¹⁾ Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be reached by reducing the speed.
2) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.

⁷⁾ When using these tubes, bucket 4255 cannot be closed with its lid.

ROTORS AND ACCESSORIES



Swing-out rotors 4176 / 4174

capacity in ml	25	3	0	1	0	2	50	2	90	650	750	1.000	1.000	175	175
Ø x L in mm	25×90	25×	:110	16:	16×80		122	62>	137	97×139	97×152	96×176	98×138	61×118	61 x 144
Cat. No.	-		-		-	512	27 ⁸⁾	-	8)	0554 ⁸⁾	0512 ⁸⁾	4239 ⁸⁾	4255	Falcon®	Nalgene®
										7)	7)	7)		7)	7)
Carrier Cat. No. 4579-A															
Cat. No.		4438		44	34		4443			-		-	4440	4430	
boring Ø x L in mm		26x72		17.5	x 53		62	x 92			98 x 138		-		-
tubes per rotor		42 28		114	1 76		6	۱ 4				6	1 4		
max. RCF ²⁾	5,8	818 <mark>5,0</mark>	26	5,886	5,094	6,113	5,320	6,113	5,411			6,294	5,501		
radius in mm	:	257 <mark>222</mark>)	260	225	270	235	270	239			278	1 243		
run-up in sec								125	125						
run-down in sec, braked								197	197						
temperature in °C1)								+4	-11						

capacity in ml	200	225	250	500
Ø x L in mm	60×130	61 x 137	60×162	96×147
Cat. No.	Nunc®	Falcon®	Cori	ning
		7)	7)	7)
Carrier Cat. No. 4579-A				
Cat. No.	4430	4440	4430	4449
boring Ø x L in mm			=	
tubes per rotor		6	1 4	
max. RCF ²⁾		6,294	5,501	
radius in mm		278	243	
run-up in sec		125	125	
run-down in sec, braked		197	197	
temperature in °C1)		+4	-11	

 $^{^{8)}}$ At temperatures of over + 40 $^{\circ}$ C and / or when not filled to capacity, bottles may warp during centrifugation.

In the following tables, the values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background.



carrier 4522-A with adapters for Falcon® tubes, blood collection tubes and bottles

Swing-out rotors 4176 / 4174

capacity in ml	4	5	6	7	12	15	25	50	100	2,6-3,4	2,7-3	4-5,5
Ø x L in mm	10x88	12x75	12×82	12×100	16x101	17×100	24×100	34×100	40×115	13×65	11×66	15x75
Cat. No.	-	0553 ²⁾	0501 ²⁾	0578 ²⁾	-	0518 ²⁾	0519 ²⁾	0521 ²⁾	0523 ²⁾	blood	collection	tubes
carrier Cat. No. 4522-A												
Cat. No.	4224	4213-93	42	13	4223	4214	4215	4216	4218	4222-93	4213-93	4214-93
boring Ø x L in mm	11 x 74	12.5 x 36	12.5	x74	16x74	17,5×74	26x74	35 x 74	41.5×74	13.2×36	12.5 x 36	17.5 x 36
tubes per rotor	378 252		288 I <mark> 192</mark>		150 100	180 120	66 44	36 I <mark>24</mark>	24 16	180 120	288 192	180 120
max. RCF ²⁾						5,818	5,003					
radius in mm						257	221					
run-up in sec						125	125					
run-down in sec, braked						197	197					
temperature in °C1)						+14	-1					

capacity in ml	4,5-5	7,5-8,5	9-10	10	1,6-5	4-7	4-7	8	8,5-10	15	50
Ø x L in mm	11x92	15×92	16×92	15×102	13x75	16×75	13×100	16x125	16×100	17×120	29x115
Cat. No.				blood co	llection/ur	ine tubes				0509	0513
							J				
carrier Cat. No. 4522-A											
Cat. No.	4213	4214	4220	4214	4222-93	4214-93	4222	4223	4214	4232	4245-A
boring Ø x L in mm	12.5x74		17.5 x 74		13.2×36	17.5x36	13.2×74	16x74	17.5×74	17×70	30×70
tubes per rotor	288 192	180 120	96 I <mark>64</mark>		180	120		150 <mark>100</mark>	180 120	138 92	48 <mark>32</mark>
max. RCF ²⁾					5,818	5,003				5,999	5,184
radius in mm		257 I <mark>221</mark>									229
run-up in sec		125 125									
run-down in sec, braked		197 197									
temperature in °C1)					+14	ı -1					

12

Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be reached by reducing the speed.
 Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote ²⁾ is 4,000.

Swing-out rotors 4176 / 4174

capacity in ml	12	25	3	0	50	250	290	600	650	750	750
Ø x L in mm	17×100	25×90	25 x	110	29x115	61 x 122	62 x 137	93×134	97 x 139	97 x 152	96×135
Cat. No.	-	-	-	-	-	5127 ⁸⁾	-	0551 ⁸⁾	0554 ⁸⁾	0512 ⁸⁾	4234-A
		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	1-								
carrier Cat. No. 4522-A											
Cat. No.	4220	4241	42	15	4249	4238		4233	4258		
boring Ø x L in mm	17.5×74	26×73	26:	ĸ74	30×96	62:	x90	94×105	97.5 x 105		
tubes per rotor	96 64	48 32	66	44	36 I <mark>24</mark>			6	١ 4		
max. RCF ²⁾		5,818	5,003		5,909 <mark>5,094</mark>	5,818	5,003		5,999	5,184	
radius in mm		257	221		261 <mark>225</mark>	257	1 221		265	1 229	
run-up in sec		125 125									
run-down in sec, braked		197 I <mark>197</mark>									
temperature in °C1)					+14	-1					

capacity in ml	250	50	500				
Ø x L in mm	60×162	96×147					
Cat. No.	Corning	Cor	ning				
carrier Cat. No. 4522-A							
Cat. No.	6322	4258 + 4449					
boring Ø x L in mm	61 x 125	97.5 x 105					
tubes per rotor	6	4					
max. RCF ²⁾	5,818 5,003	5,999	5,184				
radius in mm	257 <mark>221</mark>	265	229				
run-up in sec	125 125						
run-down in sec, braked	197 197						
temperature in °C1)	+14 -1						

Centrifugation of toxic or infectious samples:

Bucket 5052 with lid 5057 is available for use with hazardous samples. Lid 5057 provides aerosol-proof containment 9) and thus protects both the user and the environment against contamination.

Buckets 5052 and 4522-A have the same internal dimensions. This allows centrifugation of the adapters and sample tubes from the 4522-A bucket in bucket 5052 under the same conditions without the development of aerosols. Only if adapter 6322 is used will it not be possible to close bucket 5052 with lid 5057.



lid 5057

 $^{^{8)}}$ At temperatures of over + 40 $^{\circ}$ C and / or when not filled to capacity, bottles may warp during centrifugation. $^{9)}$ Tested by the TÜV in conformity with DIN EN 61010, section 2-020.

In the following tables, the values for 6-place rotor 4176 have a light turquoise background, and those for 4-place rotor 4174 have a light orange background.



rotor 4176, illustrated with carriers 4547-B and lids 5621

Swing-out rotors 4176 / 4174

capacity in ml	5	7	9	15	25	50	100	2.6-3,4	2.7-3	4-5.5	4.5-5	4.9
Ø x L in mm	12x75	12×100	14×100	17×100	24×100	34×100	44×100	13×65	11×66	15×75	11x92	13x90
Cat. No.	0553 ²⁾	0578 ²⁾	0500 ²⁾	0518 ²⁾	0519 ²⁾	0521 ²⁾	0526 ²⁾	ŀ	olood co	llection / uri	ine tubes	3
lid Cat. No. 5621												
carrier Cat. No. 4547-B					+ 0726							
Cat. No.	44	33	44	34	4438	4439	4442	4435	4433	4434	4433	4435
boring Ø x L in mm	13>	×58	17.5	x53	26×72	36×79	45x78	13.5 x 58	13×58	17.5×53	13x58	13.5 x 58
tubes per rotor	180	120	114	1 76	42 28	24 16	12 8	1261 84	180 <mark>120</mark>	114 <mark>76</mark>	180 <mark>120</mark>	126 84
max. RCF ²⁾	5,841	5,026	5,977	5,184	5,728 4,913	5,773 <mark>4,981</mark>	5,750 4,958	5,841	5,026	5,977 <mark>5,184</mark>	5,841	5,026
radius in mm	258 222 264 229 253			253 217	255 220	254 219	258	222	264 229	258	222	
run-up in sec	125 125											
run-down in sec, braked		197 I <mark>197</mark>										
temperature in °C1)						+9 1 -9)					

capacity in ml	7.5-8.5 9-10 10			1.6-5	4-7	4-7	8	8.5-10	15	50	25	;	30	10
Ø x L in mm	15x92 16x92 15x102		15 x 102	13x75	13×100	16x75	6x75 16x125 16x100		17x120	29×115	25×90	25×110		16×80
Cat. No.			blood	collectio	n / urine	tubes			0509	0513	-	-	-	-
lid Cat. No. 5621				Ī		J								
carrier Cat. No. 4547-B														
Cat. No.		4434		44	35 4434			4437	4441	4438		4434		
boring Ø x L in mm		17.5 x 53		13.5×58		17.5 x 53		17x88	30×87		26x72		17.5 x 53	
tubes per rotor		114 <mark>76</mark>		126	84		114 <mark>76</mark>		72 <mark>48</mark>	30 20		42 28		114 <mark>76</mark>
max. RCF ²⁾	5,9	977 <mark> 5,1</mark>	84	5,841	5,026	5,	977 <mark> 5,1</mark>	84	6,090	5,298	5,	909 5,	117	5,977 <mark>5,184</mark>
radius in mm	264 229 258 222					264 <mark> 229</mark>)	269	234	2	261 22	16	264 I 229	
run-up in sec		125 125												
run-down in sec, braked		197 I 197												
temperature in °C1)							+9	I -9						

Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be reached by reducing the speed.
 Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote ²⁾ is 4,000.

14

ROTORS AND ACCESSORIES



rotor 4176, illustrated with carriers 4547-B and lids 5621

Angle rotor, 6-place

 $n = 6,000 \text{ min}^{-1}$ max. RCF 5,594

₹25°

with bio-containment 9)

Swing-out rotors 4176 / 4174

Cat.	No.	4570
------	-----	------

capacity in ml	250	290	650	750	1.000	175	175	200	225	250	500		
Ø x L in mm	61 x 122	62 x 137	97 x 139	97 x 152	96×176	61 x 118	61 x 144	60×130	61 x 137	60×162	96×147		
Cat. No.	5127 ⁸⁾	- 8)	0554 ⁸⁾	0512 ⁸⁾	4239 ⁸⁾	Falcon®	Nalgene®	Nunc®	Falcon®	Cori	ning		
lid Cat. No. 5621													
carrier Cat. No. 4547-B													
Cat. No.	44	43		-		4440	4430		4440	4430	4449		
boring Ø x L in mm	62)	k 92		98×141		-							
tubes per rotor						6 4							
max. RCF ²⁾	6,203	5,411					6,384 <mark>5,592</mark>						
radius in mm	274	239	282 I <mark>247</mark>										
run-up in sec			125 I 125										
run-down in sec, braked			197 I 197										
temperature in °C1)						+9 -9	,						

capacity in ml	15	25	94	10	30	50	85	94	250
Ø x L in mm	17×100	24×100	38 x 102	16x80	26×95	29×107	38×106	38 x 102	61 x 122
Cat. No.	0518 ²⁾	0519 ²⁾	0549	-	0545	0546	0547	0538 ¹⁰⁾	5127 ⁸⁾
lid with bio-containment ⁹⁾ incl.								Ĵ	
carrier Cat. No. 4570									
Cat. No.	5646	5642	5644	5641	5642	5643	56	44	-
boring Ø x L in mm	17.6x83	26×80	38.6×88	16.6×70	26×80	29×90	38.6	5x88	61.5 x 109
tubes per rotor	42	18	6	48	18	6	(3	6
max. RCF ²⁾	5,315	5,152	4,910	5,353	5,152	4,830	4,9	910	5,594
radius in mm	132 128 122 133 128 120 122 139								
run-up in sec	64								
run-down in sec, braked	69								
temperature in °C1)					-16				

 ⁸⁾ At temperatures of over + 40 °C and / or when not filled to capacity, bottles may warp during centrifugation.
 9) Tested by the TÜV in conformity with DIN EN 61010, section 2-020.
 10) The screw cap for tube No. 0538 is available with Cat. No. 0539. Dimensions with cap 38 x 106 mm.

TECHNOLOGY	ROTO SILE	NTA 630 RS						
Floorstanding centrifuge, without roto	r coo	oled						
Power supply*)	400 V 3 ~ + N	208 – 220 V +6 / -10 % 3 ~ (+ N), with internal transformer						
Frequency	50 – 60 Hz	50 – 60 Hz						
Consumption	9,700 VA	9,000 VA						
Emission, Immunity	EN/IEC 613	26-1, class B						
Max. capacity	6 x 2,000 ml							
Max. RPM (speed)	6,000) min ⁻¹						
Max. RCF	6,5	6,520						
Running time	1 sec - 999 min : 59	1 sec - 999 min : 59 sec, ∞ continuous run						
Dimensions (HxWxD)	973 x 813	x 1,015 mm						
Weight	approx	. 355 kg						
Refrigeration								
Temperature control, infinitely variable	from - 20	from - 20 to + 40 °C						
Cat. No.	5005	5005-08						

^{*)} Other voltages on request.









Hettich centrifuges comply with all relevant EU standards in effect and conform to the European level of quality and safety for medical devices. Evidence is provided by national and international test marks such as IEC 61010 or the CE conformity. The ISO 9001, ISO 13485 and ISO 14001 certificates accredited to the company bear witness to the extreme care and responsibility Hettich puts into the manufacturing of centrifuges and

Our certification as an "Authorised Economic Operator" enables accelerated customs clearance.

Andreas Hettich GmbH & Co. KG

Föhrenstr. 12 D-78532 Tuttlingen Germany www.hettichlab.com info@hettichlab.com service@hettichlab.com Phone +49 (0)7461/705-0

+49 (0)7461/705 -1122

National Sales: -1200 -1201 **International Sales: National Service:** -1202 International Service: -12 03

